

# User Manual Labcare V 3.0.1

---

## PHLEBOTOMY MODULE

**P S Solutions**

P.O. BOX 31291, AL-JAZEERA, AL-HAMRA, RAK, UAE.

## Table of Contents

1. PHLEBOTOMY .....	2
1.1 Login Screen .....	2
1.2 Sample Collection.....	3

## 1. PHLEBOTOMY

### 1.1 Login Screen



- After double clicking on Labcare Short cut available on user desktop user will see this screen as Login Screen of Labcare application.
- From center option user has to select main center, according to which users will be loaded for particular center.
- This is the login screen of Labcare software. Only authorized users can log on to the system.
- The user has to select the user name and enter the password then click on "OK" to proceed.
- If the user name or password is wrong then the user will be prompted to enter the correct username or password. If it is an authorized user then the registration screen will be opened.
- On clicking the cancel button the application ends.

## 1.2 Sample Collection

**Sample Collection**

Date: >= 16/12/2016 LabNo: Ref By: ALL Sample Status: Registered Category: < All > Stat: ☐

**Total Patient(s) [3]**

LabNo	Visit Date	IPOP No	Title	Patient Name
190000027	17/12/2016	500000007		Mus Test 01 [M35 YEARS]
190000026	17/12/2016	500000007		Mus Test 01 [M35 YEARS]
190000025	17/12/2016	500000007		Mus Test 01 [M35 YEARS]

**Comment**

Sample Type	Tubes / Containers	Collection Comment
Aerobic Blood		
Amniotic Fluid		
Anaerobic Blood		

Test Name	Category	Profiles	Stat	Processing Centre	Outsic	Select	Suffix
<input checked="" type="checkbox"/> Food Allergens Panel	ALLERGY	Package	N	ANH AIRPORT	N	<input checked="" type="checkbox"/>	Aerobic Blo
<input checked="" type="checkbox"/> Inhalation Allergens Panel	ALLERGY	Package	N	ANH AIRPORT	N	<input checked="" type="checkbox"/>	Amniotic Fl
<input checked="" type="checkbox"/> Pediatric Allergens Panel	ALLERGY	Package	N	ANH AIRPORT	N	<input checked="" type="checkbox"/>	Anaerobic E
<input checked="" type="checkbox"/> Albumin; Random urine	BIOCHEMISTRY	Package	N	ANH AIRPORT	N	<input checked="" type="checkbox"/>	Any Sample
<input checked="" type="checkbox"/> Post Prandial Glucose	BIOCHEMISTRY	Package	N	ALNOOR_KH	Y	<input checked="" type="checkbox"/>	Ascitic Body
<input checked="" type="checkbox"/> Growth Hormone Suppression	BIOCHEMISTRY	Package	N	ALNOOR_KH	N	<input checked="" type="checkbox"/>	Aspiration f
<input checked="" type="checkbox"/> Growth Hormone Stimulation	BIOCHEMISTRY	Package	N	ALNOOR_KH	N	<input checked="" type="checkbox"/>	BAL
<input checked="" type="checkbox"/> Insulin Tolerance For Growth Hormone	BIOCHEMISTRY	Package	N	ALNOOR_KH	N	<input checked="" type="checkbox"/>	Bile Body Fl
						<input checked="" type="checkbox"/>	Blood

☐ Select All ☐ Stat Sample

Buttons: Print, Reject, Close

- User can select "Sample Collection" by going to "Sample" in Labcare menu.
- This is Labcare sample collection window, once the patient is registered, it will be available in sample collection window. Without sample collection, it will not allow the tests for sample receiving (source).
- User can enter no. of Sample Tubes/Container collected from the patient and collection comment under the comment section.
- User can collect sample by Clicking in "Print" button, it will generate barcodes according to sample type and Quantity selected and this will change sample status of the patient from "Registered" to "Sample Collected".
- User can reject the test by using "Reject" button. This test will be known as "Rejected" test.
- To reprint barcode, select sample status as "Sample Collected" and then click on Reprint button.